






## PELVIC ELECTROSTIMULATORS

						
	<b>BIOLITO</b>	<b>MYOLITO</b>	<b>CONTINENCE</b>	<b>PELVIFINE</b>	<b>PELVITONE</b>	<b>EVOSTIM-UG</b>
n. channels	2	2	2	2	2	2
Stim currents	continuous work / rest (WR) modulated - burst	continuous work / rest (WR) modulated - burst	continuous work / rest (WR)	continuous work / rest (WR)	continuous work / rest (WR) modulated	continuous work / rest (WR)
Tens currents	han	han	yes	yes	yes	yes
free program	1	1	3+1 (cont)	3	3	10
n ° sessions pr. free	1/5	1/5	1	1/5	1/5	1/10
accessories	electrodes 4x40x40 - 4x40x80	electrodes 4x40x40 - 4x40x80	No	no	no	Vaginal probe MINIMA
N° programs	5	5 + 7	9	11	11	5x10
Urgency incontinence	1	1	1	1	1	1
Stress / Stress incontinence	2	2	2	2	2	1
Faecal incontinence	1	1				1
Mixed incontinence	1	1	1	1	1	1
Pain reliever			2	3	3	
Frequency / Urgency			1	1	1	
Loss of sensation			1	1	1	1
Pelvic floor (toning)				1	1	
Relaxation			1	1	1	
Endurance-Resistance		4				
Nerve stimulation		3				

## BIOFEEDBACK AND/OR ELECTROSTIMULATORS

				
	<b>SIMPLEX</b>	<b>SIMPLEX PRO</b>	<b>EVOSTIM P</b>	<b>EVOSTIM E</b>
n. channels	1	2	2	2
currents			Continue - work / rest (W / R) - ETS	Continue - work / rest (W / R) - ETS
Biofeedback	electromyography	pressure electromyography	pressure	electromyography
screen	black and white	black and white	black and white + colour backlight	black and white + colour backlight
software	yes	yes		
Incontinence programs			5	5
Accessories			Vaginal pressure probe V2STFW	Vp. PERISPHERA- O - 4 electrodes 40x40